



Case Study

Transforming Biologics Manufacturing Through Smart Upstream Process Intensification



In every molecule is the possibility for better health.

Overview

The rapidly evolving biologics industry is driven by the need for higher productivity, cost efficiency, and reduced facility footprints. Traditional fed-batch processes have reliably produced therapeutic proteins but face limitations in meeting growing market demands. Commonly explored intensification strategies include N-1 perfusion, full perfusion at production scale, and concentrated fed-batch processes. While each offers advantages, they also come with increased complexity, capital investment, or marginal gains.

In response, we developed an integrated upstream intensification strategy that combined high cell density cell banks, optimized seed train perfusion, and a fortified concentrated fed-batch process. This approach synergized the strengths of various intensification methods while minimizing their drawbacks, delivering significant productivity improvements, and maintaining product quality.

Challenge

Industry-standard approaches to intensification typically involve trade-offs. Traditional fed-batch processes offer moderate titers (~4–6 g/L) and monthly outputs around 40–50 kg using 4 x 2kL reactors. Full-scale perfusion production boosts productivity further but demands substantial capital investment and staffing resources. N-1 perfusion can shorten seed train duration and modestly increase titer, but it introduces a slightly high media consumption and increases operational complexity. Concentrated fed-batch processes improve titers but require more media than N-1 perfusion (Figure 1). Each approach offered incremental gains—but also came with trade-offs in cost, complexity, or scale-up risk. So, the choice of intensification is very important to get a meaningful outcome.




Figure 1: Comparison of upstream process intensification scenarios, illustrating traditional fed-batch bioprocessing alongside advanced approaches such as N-1 intensification and perfusion-enabled strategies.

The challenge was to design an upstream intensification platform that achieves the productivity benefits of advanced perfusion systems without excessive complexity or capital expenditure, while ensuring appropriate product quality.

Aragen's Approach

To accelerate the overall bioprocess timeline while maintaining cell quality and productivity, a streamlined and intensified seed train and production strategy was implemented. The process (Figure 2) began with the CSPR screening at shake flasks to identify best condition to achieve a higher viable cell density at shake flask, establishment of high cell density cell banks (30–70 million cells/mL) to provide a robust foundation for seed expansion. Subsequently, N-1 perfusion was performed using a perfusion device, enabling rapid expansion of viable cell numbers. Upon reaching the desired cell density, N stage production-scale bioreactor was directly inoculated from the N-1 stage, bypassing intermediate bioreactor steps to streamline the workflow.

Process Intensification Workflow




Figure 2: Process intensification workflow.

To further optimize scale-up, traditional intermediate-scale bioreactors could be replaced with perfusion-enabled Wave bioreactors at the seed stage, simplifying process logistics and reduce Cost of Goods Sold (COGS).

At Aragen, for process intensification, a concentrated fed-batch approach has been chosen to reduce the operation complexity and betterment of COGS. In this process a 50 kDa hollow fibre membrane has been used to retain monoclonal antibodies and allowing only waste metabolites to be removed from the system while ensueing the continuous addition of fresh nutrients.

Fresh media addition rate was adjusted in such a way that it helps to increases the cell mass. The primary objective was to get a higher a viable cell density (VCD) of > 100 million cells/mL in a short period. Once the viable cell density reaches >100 million cells/mL, 75% culture is harvested and purified up to Protein A, followed by LPT and intermediate depth filtration steps. At this stage, spent media analysis should also be performed to understand the key media component depletion trends of key media components. At Aragen, a robust spent media analysis method has been developed, which can identify 85 components divided into seven categories: amino acid, acids, vitamins, amino acid derivatives, nucleotides, nucleosides, choline-related components, and others.

A fresh media top-up will then be performed with the remaining 25% culture to achieve the final VCD in a range of 25-30 million cells/mL. The perfusion rate will be adjusted to maintain a VCD between 50-60 million cells/mL, providing greater operational flexibility using a single-use bioreactor. At this stage, spent medium analysis data should be leveraged to fortify the media with required amino acids, reducing media consumption while maintaining optimum nutrient levels.



Outcomes

The intensification strategy combined high cell density cell banks, N-1 perfusion, concentrated fed-batch, and targeted media fortification to enhance productivity and efficiency without compromising product quality or increasing complexity.

- **Cell-Specific Perfusion Rate (CSPR) Screening:** From the CSPR screening at shake flask scale, it is observed that 20-30 (pL/cell/day) CSPR is the best fit for the process (Figure 3). This understanding is helpful to determine the optimal perfusion rate at N-1 and N stage.




Figure 3: Cell-Specific Perfusion Rate (CSPR) screening at shake flask scale. CSPR between 20-30 (pL/cell/day) was ideal for this process.





- **High Seed Density Cell Bank (HSDCB):** Simultaneously, a high-density cell bank is prepared in regular cryovials at a maximum concentration of 70 million cells/mL. The thaw viability and the post-thaw doubling time of the HSDCB have been evaluated, and results indicate excellent recovery and robust cell growth in subsequent passages (**Figure 4**).




Figure 4: Doubling time and thaw viability across High Seed Density (HSD) conditions. Bars show doubling time (hrs); orange diamonds indicate thaw viability (%).

- **Seed Train and Perfusion:** HSDCB reduced seed train time by 30–40%. By combining the HSDCB and N-1 perfusion, at N-1 stage reached ~40 million cells/mL by Day 5 (Figure 5).




Figure 5: Vessel Volume Per Day (VVD) and Viable Cell Concentration (VCC) for N-1 perfusion process.

- Concentrated Fed-Batch Performance:** At first with higher VVD, cell densities up to 124 million cells/mL and titer of 12 g/L by Day 11 (Figure 6A) were achieved with stable productivity (qP ~20 PCD) (Figure 6B).




Figure 6: Volumetric Exchange Rate (VVD) and Viable Cell Concentration (VCC) for concentrated fed-batch process (A) and specific productivity (qp in pg/cell/day) (B).

- Media Fortification:** Once the VCD reached to 124 million cells/mL, 75% of the broth was harvested. The remaining 25% was reinitiated with fresh media, starting at an initial viable cell density of 28 million cells/mL. The perfusion media was fortified to maintain the required amino acids to support the cell growth. Targeted nutrient supplementation sustained an average cell density of 45 million cells/mL and boosted titers to 15 g/L by Day 10 at low media usage (VVD 0.7), while maintaining consistent productivity (qP ~40 PCD) (Figure 7).




Figure 7: Volumetric Exchange Rate (VVD), Viable Cell Concentration (VCC) and titre (in g/L) for concentrated fed-batch with perfusion media fortification process (A) and specific productivity (qp in pg/cell/day) (B).



- N-Stage Intensification:** Cumulative titres reached 27 g/L (>3× traditional fed-batch) (Figure 8A); daily productivity increased to 1.3 g/L/day. Using this approach at 4x2kL bioreactors monthly output can be achieved ~140 kg (Figure 8B).




Figure 8: N-stage intensification with intermediate harvests resulted in a threefold increase in titer, achieving a cumulative productivity of 27 g/L by day 21 **(A)** Bar graph comparing titers (g/L/day) of three processes-traditional fed batch (nFB), High Seed Density (HSD) and Intensified Fed Batch (iFB) **(B)**.

- Product Quality:** Glycan profiles and post-translational modifications remained consistent with biosimilar standards, ensuring regulatory compliance (Figure 9).




Figure 9: Glycan profile comparison across biosimilars and intensified fed-batch samples.

Blue bars represent Biosimilars 1–3, and red bars represent Fed-batch samples (Days 9, 11, 19, and 21), compared across GOF, Total Galactosylation, and Terminal Glc-Nac categories. The glycan profiles show strong consistency, with aligned GOF levels, acceptable variation in galactosylation, and matched Glc-Nac levels, confirming structural and functional comparability.

This approach delivered a threefold increase in titer and production capacity while maintaining cell viability and product quality, all without the higher costs of full perfusion systems.

Conclusion

By integrating high-density cell banking, optimized perfusion seed trains, and fortified concentrated fed-batch culture, our approach takes upstream process intensification beyond traditional methods. This strategy allows biologics manufacturers to achieve much higher product titres and monthly output, all while keeping complexity and costs manageable.

Our results demonstrate a practical and scalable pathway for next-generation biologics manufacturing—one that balances productivity, quality, and cost efficiency. Importantly, this approach can be implemented without large-scale perfusion facilities, helping manufacturers meet growing global demand, use resources more efficiently, and maintain regulatory compliance.

[Click Here to Explore More Case Studies](#)



www.aragen.com

Aragen is a trusted R&D and manufacturing partner to the global life sciences industry. From concept to commercial, we transform your ideas into solutions for better health. We have more than 20 years of experience providing a range of contract research, development and manufacturing services across the drug development continuum, with a focus on early-stage discovery and development of new molecular entities (NMEs). Whether large pharma or biotech, an agrochemical or animal health company, we provide you global resources and proven capabilities at every stage of the biopharma lifecycle, in small and large molecules. Our ability to offer end-to-end solutions or support standalone programs is underpinned by an innovation mindset, enabling technologies, and a partnership approach to every engagement.

At Aragen, we recognize your work is vital, urgent and impacts lives. Our purpose, 'In every molecule is the possibility for better health' motivates us to drive the success of your programs, so that we can together transform hope into health for millions of people around the world.

For more information, contact info@aragen.com

© 2025 Aragen Life Sciences, All Rights Reserved.

Stay Connected

